

## AUTOCLAVABLE CUFF ARM CONICAL 65X12CM

Ref. GHP552804

# Hpm™


CE



 30/08/2024 09:57:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## AUTOCLAVABLE CUFF ARM CONICAL 65X12CM



Ref. GHP552804

**Operating room equipment: double pocket autoclavable cuff for two-way electronic pneumatic tourniquet. Sold individually.**

### Description

Nylon cuff with single silicone bladder suitable for autoclaving. With male Luer lock connector and cap for cleaning and sterilisation.

### Features

CE Marking	CE	Medical device	No
Bar code	3700446033682	HS Code	9018 90 2090
Warranty (months)	12	Country	Sweden
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.32
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 09:57:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)