

AUTOCLAV.DOUBLE AD CONICAL ARM CUFF(47X16CM) W/CONNECTOR F/EGP200

Hpm™


Ref. GHP552813



 30/08/2024 10:00:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

AUTOCLAV.DOUBLE AD CONICAL ARM CUFF(47X16CM) W/CONNECTOR F/EGP200



Ref. **GHP552813**

Operating room equipment: double pocket autoclavable cuff for two-way electronic pneumatic tourniquet. Sold individually.

Description

Nylon cuff with double silicone bladder suitable for autoclaving. With male Luer lock connector and cap for cleaning and sterilization.

Features

CE Marking	CE	Medical device	No
Bar code	3700446033774	HS Code	9018 90 2090
Warranty (months)	12	Country	Sweden
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.615
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 10:00:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr