

AUTOCLAVABLE DOUBLE CUFF SILICON BLADDER SMALL CHILD 46X10CM

Ref. GHP552817

Hpm™


CE



 30/08/2024 09:58:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

AUTOCLAVABLE DOUBLE CUFF SILICON BLADDER SMALL CHILD 46X10CM



Ref. **GHP552817**

Operating room equipment: double pocket autoclavable cuff for two-way electronic pneumatic tourniquet. Sold individually.

Description

Nylon cuffs with double autoclavable silicone pocket.
With male Luer Lock tip and cap for cleaning and sterilization.
For orthopedic surgery and traumatology.
No risk of latex allergy.
Velcro and nylon clamping strap to increase safety.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446033811	HS Code	9018 90 2090
Country	Sweden		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.196
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 09:58:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr