

DISPOSABLE SINGLE CUFF 109X13CM GREEN (CARTON OF 10 PIECES)

Ref. GHP552898

Hpm™

CE
0123

DISPOSITIF MEDICAL
CLASSE
Is
MEDICAL DEVICE

GARANTIE
1
YEAR
GUARANTEE



📅 30/08/2024 09:59:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

DISPOSABLE SINGLE CUFF 109X13CM GREEN (CARTON OF 10 PIECES)



Ref. GHP552898

HPM armbands/bib shorts give you better access to the operating area thanks to a 2 m long TPU hose. The cuffs are made of high-quality material, which makes them soft and easy to apply, without the risk of latex allergy. Wide range of sizes with color coded for easy identification.

Description

ETO Sterilized. 2 m long hose made of TPU. Straight and Conical. No Latex. Wide range of sizes. Color coded. Customized protection sleeve included. CPC male as standard. Soft and easy to apply

Features

| | | | |
|------------------------------------|---------------|---------------------------|-----------------------------|
| CE Marking | CE0123 | Medical device | Class Is (sterile function) |
| Bar code | 7350054600279 | HS Code | 9018 90 2090 |
| Warranty (months) | 12 | Country | Sweden |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 26 | Packed unit / Height (cm) | 20 |
| Packed unit / Length (cm) | 64 | Packed unit / Weight | 3.858 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

30/08/2024 09:59:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr