

DISPOSABLE SPECULUM DIA. 2.5 MM, BOX OF 250, BLACK, CHILD, SPENGLER



Ref. L20747



30/08/2024 21:34:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

DISPOSABLE SPECULUM DIA. 2.5 MM, BOX OF 250, BLACK, CHILD, SPENGLER



Ref. L20747

Description

Box of 250.
For children, single use,
Black,
Ø 2.5 mm
Universal, compatible with all brands
Box with dispensing spout for easy access to speculums

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446031268	HS Code	90330000
Country	Sweden		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	9	Packed unit / Height (cm)	11
Packed unit / Length (cm)	11	Packed unit / Weight	0.14
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 21:34:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr