


MINI BAG EIFFEL, DE BOISSY

Ref. MDEIFF02

De Boissy
— DEPUIS 1982



 30/08/2024 22:15:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

MINI BAG EIFFEL, DE BOISSY

Ref. MDEIFF02

De Boissy
— DEPUIS 1982

The mini version of the EIFFEL Medical Bag

Description

Chic and simple at the same time, the EIFFEL mini bag recalls the shape of the traditional leather bag of doctors.

Thanks to its bidirectional opening and multiple pockets/compartments, it makes it easy to store and reach all the care equipment and diagnostic instruments that healthcare professionals need.

Of a mottled gray, the polyester fabric is resistant, durable, water-repellent and therefore washable, it is accompanied by a reflective strip for safety, and YKK® slides for their quality.

- A main space that includes 6 pockets and 1 removable velcro divider. The large bidirectional zip opening gives very easy access to diagnostic instruments and care equipment
- A ventral pocket (31x15cm) opens completely with on one side 2 zipped storage nets and 2 flat pockets, and on the other side several elastic loops that allow to store small care supplies: pliers, scissors, pen lamp, thermometer, bandages, pens, etc...
- There are also four small external zipped pockets and another on the side with elastic to slide the bottle of hydroalcoholic gel that remains easily accessible at any time
- Padded and solid handles, a padded and adjustable shoulder strap - A name tag to identify the bag
- A non-slip bottom to reduce friction and wear
- Dimensions w x h x d: 32 x 28 x 20 cm- Weight: 0.95 kg
- Maximum recommended content: 4 kg / 15 liters
- Sold empty

 30/08/2024 22:15:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France


 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

MINI BAG EIFFEL, DE BOISSY

Ref. MDEIFF02



Features

CE Marking	NONE	Medical device	No
Bar code	3661040218946	HS Code	4202220000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	28	Bare Unit / Width (cm)	28
Packed unit / Height (cm)	8	Bare Unit / Height (cm)	20
Packed unit / Length (cm)	42	Bare Unit / Length (cm)	32
Packed unit / Weight	nc	Bare Unit / Weight (kg)	0.94
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 22:15:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr