

## CONVERTIBLE BACK PACK MADRAGUE, DE BOISSY, WITH BLACK LEATHER POCKET

*De Boissy*  
— DEPUIS 1982

Ref. MDMADR03



 30/08/2024 22:33:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

# CONVERTIBLE BACK PACK MADRAGUE, DE BOISSY, WITH BLACK LEATHER POCKET

*De Boissy*  
— DEPUIS 1982

Ref. **MDMADR03**

**It is worn in a backpack, hand or shoulder; A real Swiss army knife!**

## Description

Eco-responsible, customizable and "super organized" for a trendy look with all the practical and modern options.

All fabrics, exterior and interior, certified Global Recycled Standard (GRS), have been made from recycled material and remain washable, resistant and water-repellent.

Madrague : it is the pink-red of the laurels of the Mediterranean, a symbolic place, refuge beach of a female icon of French cinema of the "sixties".  
Madrague bags have been designed to be "ultra practical", to be the best ally of health professionals during their medical visits at home or in community.  
- Its main space

is modular and contains 4 zipped compartments, extractable (fixed with velcro): 3 transparencies to quickly identify the contents, and 1 isothermal, supplied with its reusable cold gel pocket, to slip up to 20 bulbs

-An extension zip makes it possible to enlarge this main compartment to store more equipment

- A second padded space can store computer / tablet and documents

- This model also has the large ventral pocket (45x21 cm) which opens completely with on one side 2 zipped storage nets and 2 flat pockets, and on the other side several elastic loops that allow to store small care supplies: pliers, scissors, pen lamp, thermometer, bandages, pens, etc ...

-Two external zipped pockets, plus one with elastic dedicated to the bottle of hydroalcoholic gel, a pouch, detachable-customizable grained leather appearance (printed fancy pouch sold separately)

- An external USB port with integrated charging cable allows you to connect a backup battery and never run out of phone, or forget its charger at a patient

- The shoulder straps, handles and shoulder strap are padded, solid and concealable to adopt the most suitable wear

-Dimensions w x h x d: 30 x 45 x 15/19 cm- Weight: 2.2 kg

-  
Maximum recommended content: 6 kg / 19-25 liters

The De Boissy bag is put back in its dust bag to protect it when stored. Grey with patterns, this accessory is designed from recycled materials with a cord tightening.

 30/08/2024 22:33:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## CONVERTIBLE BACK PACK MADRAGUE, DE BOISSY, WITH BLACK LEATHER POCKET



Ref. **MDMADR03**

### Features

CE Marking	NONE	Medical device	No
Bar code	3661040219011	HS Code	4202129000
Warranty (months)	12	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	50	Bare Unit / Width (cm)	15
Packed unit / Height (cm)	35	Bare Unit / Height (cm)	45
Packed unit / Length (cm)	50	Bare Unit / Length (cm)	30
Packed unit / Weight	1.86	Bare Unit / Weight (kg)	1.86
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 22:33:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)