

## DISPLAY OF 12 THERMOMETRES TEMPO 10 BLUE



Ref. SP1623P



30/08/2024 21:25:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## DISPLAY OF 12 THERMOMETRES TEMPO 10 BLUE



Ref. **SP1623P**

### Description

Display of 12 TEMPO 10 blue thermometers. The TEMPO 10 thermometer is a thermometer showing body temperature in 10 seconds\*

- Fever alarm\*
- Automatic shutdown
- Water resistant mouthpiece\*
- Last measurement memory\*
- Measurement in 10 seconds
- Comes with a storage case and button cell battery (1.5V)

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446027513	HS Code	90251920
Warranty (months)	24	Country	China

#### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.65

#### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

#### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

30/08/2024 21:25:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)