

## DISPLAY OF 12 THERMOMETRES TEMPO 10 FLEX GREY



Ref. SP1648P



30/08/2024 21:26:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## DISPLAY OF 12 THERMOMETRES TEMPO 10 FLEX GREY



Ref. **SP1648P**

### Description

Display of 12 thermometers TEMPO 10 Flex gray The thermometer TEMPO 10

FLEX is a thermometer indicating the body temperature in 10 seconds and having a flexible

mouthpiece• Fever alarm•

Automatic shutdown•

Water-resistant and flexible tip•

Last measurement memory

• Measurement in 10 seconds•

2 colors available green and gray• Comes with a storage case and button battery (1.5V)

- Comes with: Storage case, button battery (1.5V).

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446041625	HS Code	90251920
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	16
Packed unit / Length (cm)	11	Packed unit / Weight	0.66
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 21:26:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)