

## THERMOMETER TEMPO 10 FLEX GREY

Ref. SP1648



 30/08/2024 21:27:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## THERMOMETER TEMPO 10 FLEX GREY



Ref. **SP1648**

### Description

The TEMPO 10 FLEX thermometer is a thermometer indicating body temperature in 10 seconds and has a flexible mouthpiece•

- Fever alarm•
- Automatic shutdown•
- Water resistant and flexible mouthpiece•
- Last measurement memory•
- Measurement in 10 seconds•
- 2 colors available green and gray
- Comes with a storage case and button battery (1.5V)

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446041618	HS Code	9025192090
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	6	Packed unit / Height (cm)	2
Packed unit / Length (cm)	15	Packed unit / Weight	0.052
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 21:27:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)