

## SPO2 SENSOR WRAP FOR MASTER PALM 2



Ref. 222253



30/08/2024 16:46:15 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## SPO2 SENSOR WRAP FOR MASTER PALM 2



Ref. 222253

### Description

SpO2 wrap sensor for oximeter MASTER PALM 2

### Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446030230	HS Code	90181990
Warranty (months)	12	Country	Uruguay
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.085
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	1	Max temperature	50
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 16:46:15 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)