

## SILICON PROTECT COVER FOR MASTER PALM 2/3 (2 PCS)



Ref. 223004



CE  
0123



 30/08/2024 16:54:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## SILICON PROTECT COVER FOR MASTER PALM 2/3 (2 PCS)



Ref. 223004

### Description

Silicone protection for MASTER PALM 2 and MASTER PALM 3 oximeters.  
Set of 2 protections, one yellow and one transparent.

### Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	3700446025434	HS Code	90330000
Country	China		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.074
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 16:54:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)