

CARDIOMATE PC ECG SPENGLER

Ref. 260039



30/08/2024 21:52:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CARDIOMATE PC ECG SPENGLER



Ref. 260039

A SPENGLER PC-ECG with XXL performance in the palm of your hand

Description

Digital ECG, the Cardiomate PC allows the real-time display of a 12-lead ECG on the screen of a computer and its printing on an A4 printer. Compatible with a Windows system, the Cardiomate PC connects via a simple USB cable. With the same performance as the Cardiomate range and with its compact and lightweight case, the Cardiomate PC is a state-of-the-art measuring instrument in the office or on the go. The software allows the visualization, recording, analysis and storage of ECGs in a clear and intuitive way.

VERY HIGH ACCURACY OF MEASUREMENTS- Visualization of the 12 leads in real time on PC
 -Sampling speed 32,000Hz
 -24-bit converter
 -Automatic interpretation and calculation of intervals
 -Pacemaker detection

EASY TO USE
 -Connection to the PC by USB
 -Self-powered- Simple and intuitive software
 -Software compatible with other ECG Cardiomate

COMPACT AND LEGER
 -Dimensions: 68x126x24mm
 - Weight: 104g

Features

| | | | |
|------------------------------------|---------------|-------------------------|------------|
| CE Marking | CE0123 | Medical device | Class IIa |
| Bar code | 3700446042059 | HS Code | 9018110000 |
| Warranty (months) | 24 | Country | China |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Bare Unit / Width (cm) | 12.6 |
| Packed unit / Height (cm) | nc | Bare Unit / Height (cm) | 2.4 |
| Packed unit / Length (cm) | nc | Bare Unit / Length (cm) | 6.7 |
| Packed unit / Weight | nc | Bare Unit / Weight (kg) | 0.13 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

 30/08/2024 21:52:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
 13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr