

BLACK BEZEL FOR SPENGLER ANEROID SPHYGMOMANOMETERS LIAN NANO AND VAQUEZ-LAUBRY NANO




Ref. 310163



 30/08/2024 17:21:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

BLACK BEZEL FOR SPENGLER ANEROID SPHYGMOMANOMETERS LIAN NANO AND VAQUEZ-LAUBRY NANO



Ref. 310163

Description

Spare part SPENGLER, Black bezel for SPENGLER aneroid sphygmomanometers Lian Nano and Vaquez-Laubry Nano

Features

CE Marking	NONE	Medical device	No
Bar code	3700446028626	HS Code	9018901000
Warranty (months)	60	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 17:21:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr