

GLASS FOR SPENGLER VAQUEZ LAUBRY COLOUR/HOPITAL (LARGE DIAL) (X12)




Ref. 310220



 30/08/2024 17:22:10 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

GLASS FOR SPENGLER VAQUEZ LAUBRY COLOUR/HOPITAL (LARGE DIAL) (X12)



Ref. 310220

Description

Kit of 12 replacement lenses for blood pressure monitors cuff VAQUEZ-LAUBRY HOSPITAL

Features

CE Marking	NONE	Medical device	No
Bar code	3700446024901	HS Code	901890100000L
Warranty (months)	60	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	2	Packed unit / Height (cm)	2
Packed unit / Length (cm)	2	Packed unit / Weight	1.26
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 17:22:10 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr