

## CHROME-PLATED LAUBRY STETHOSCOPE / BLACK



Ref. 507300



CE



30/08/2024 21:03:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CHROME-PLATED LAUBRY STETHOSCOPE / BLACK



Ref. 507300

=====

### Description

The famous cardiophone. Single head chestpiece (107 g) made of heavy brass. Elegant and streamlined for easy handling. High resolution, fiberglass, acoustic diaphragm, angled binaural and screw-on flexible tips. Ideal for monitoring blood pressure and general auscultation.

### Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009151	HS Code	90181990
Warranty (months)	60		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	32	Packed unit / Weight	0.393
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Colors	Black	Packed unit / Quantity	1

30/08/2024 21:03:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)