

## DIGITAL OTOSCOPE FIREFLY DE550

Ref. DFIDE550



30/08/2024 21:51:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## DIGITAL OTOSCOPE FIREFLY DE550



Ref. DFIDE550

Connected otoscope for telemedicine use

### Description

High precision digital handheld otoscope video WIRELESS; makes very precise and clear observations of the ear.

- Resolution of 720x480 real pixels without processing with interpolation techniques.
- Provides instant images and videos: DE500 uses a state-of-the-art image processing engine that virtually eliminates visual lag, allowing for an intuitive user experience.
- Magnification up to 50x optically and 150x digitally for accurate analysis.
- Equipped with multi-layered glass lenses, rather than plastic, thus offering superb depth of field and unparalleled image sharpness.
- High-quality ultra-bright LEDs with adjustable brightness are installed.
- Easily integrates with digital medical records (DME)

Sensor Resolution: 720 x 480- Image - Video:

Format: YUY2 / Frame rate: 30 FPS / Format: JPEG / BMP / AVI- Image

properties / Video: Color: Saturation, Hue / Exposure: Brightness, Contrast / Image: Image sharpness, Gamma- Fully adjustable brightness: 4 ultrabright LEDs / Fully adjustable brightness

- Accessories: Speculums 3/4/5 mm included

- Windows 10, 8, 7, and Mac (Intel) OS-X 10.4+ operating systems

### Features

CE Marking	CE	Medical device	Class I
Bar code	8512890070190	HS Code	90189042
Warranty (months)	12		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	20	Packed unit / Height (cm)	5
Packed unit / Length (cm)	15	Packed unit / Weight	0.302
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 21:51:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)