

SUNTECH CUFF SINGLE PATIENT SCREW TIGHT 38-50CM

SUNTECH

Ref. DSU231059



30/08/2024 14:58:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SUNTECH CUFF SINGLE PATIENT SCREW TIGHT 38-50CM



Ref. **DSU231059**

Description

Set of 20 Suntech "SOFT" single-patient armbands, brown.
Bib shorts size 38-50cm, screw connection.


Features

CE Marking	CE	Medical device	Class I
Bar code	3661041009772	HS Code	9018905040
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 14:58:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr