

SET OF 5 SUNTECH MONOBLOC ARMBANDS REUSABLE SMALL ADULT (17-25CM) - BAYONET METAL CONNECTOR

SUNTECH

Ref. DSU231112

CE



 30/08/2024 14:49:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SET OF 5 SUNTECH MONOBLOC ARMBANDS REUSABLE SMALL ADULT (17-25CM) - BAYONET METAL CONNECTOR



Ref. [DSU231112](#)

Description

Reusable nylon cuff sold in boxes of 5 units.
Single tubing, bayonette connection 17-25cm.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661041008805	HS Code	9018905040
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	18	Packed unit / Height (cm)	13
Packed unit / Length (cm)	9	Packed unit / Weight	0.4
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 14:49:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr