

OPD CUFF ADULT 23-33CM METAL MALE BAYONET (BOX OF 5)

SUNTECH

Ref. DSU231113



CE



 30/08/2024 14:50:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

OPD CUFF ADULT 23-33CM METAL MALE BAYONET (BOX OF 5)

SUNTECH

Ref. **DSU231113**

Description

Reusable nylon cuff sold in boxes of 5 units.
Single tubing, bayonette connection 23-33cm.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661041008812	HS Code	9018905040
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	18
Packed unit / Length (cm)	9	Packed unit / Weight	0.45
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 14:50:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr