



EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

EMERGENCY'S™

Ref. MEEMER02_RC

GARANTIE
2 ANS
2 YEARS
GUARANTEE



 30/08/2024 08:01:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

EMERGENCY'S™

Ref. **MEEMER02_RC**

Practical and lightweight first aid bag with two adjustable internal dividers and a zippered pocket in the flap to facilitate the organization of sanitary equipment

Description

- It also has three outer pockets; the one at the front is divided with rubberized elastic bands
- High resistance zip closures
- Reflective strips- Washable case
- Reinforced handles and shoulder straps
- Dimensions 44 x 25 x 27 cm, capacity 22 L
- Material 600D Polyester
- Weight 1.29 kg
- Available in waterproof version
- Sold empty

Features

CE Marking	NONE	Medical device	No
Bar code	8436564542670	HS Code	420212
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	24	Bare Unit / Width (cm)	25
Packed unit / Height (cm)	22	Bare Unit / Height (cm)	27
Packed unit / Length (cm)	30	Bare Unit / Length (cm)	44
Packed unit / Weight	1.3	Bare Unit / Weight (kg)	1.29
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Pink	Packed unit / Quantity	1

 30/08/2024 08:01:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr