

## NON-CONTACT THERMOMETER TEMPO PRO SPENGLER



Ref. SP3110



CE  
0123

DISPOSITIF MEDICAL  
CLASSE  
**IIa**  
MEDICAL DEVICE

GARANTIE  
**2**  
2 ANS  
GUARANTEE



30/08/2024 21:39:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## NON-CONTACT THERMOMETER TEMPO PRO SPENGLER



Ref. SP3110

**Infrared (IR), non-contact, clinical thermometer, digital, handheld, range 34 # 43 °C**

### Description

At a distance of a few centimetres from the temple, its display instantly shows the temperature.

Without contact to avoid cross-contamination, Tempo Pro is easy to disinfect and always stays clean.

#### TECHNICAL CHARACTERISTICS

- Measurement : Infrared
- Distance de mesure : - 5 cm
- Unit of measurement : Celsius or Fahrenheit
- Measuring range : 32,0°C – 43,0°C / 89,6°F - 109,4°F
- Resolution : 0,1°C / F°
- Forehead accuracy : 0,2°C / 0,4F°
- Memories : 30 memories
- Measurement time: Less than 1 second
- Screen : Backlit LCD screen
- Auto power off
- 2 x 1,5V AAA batteries (included)

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446041823	HS Code	90251920
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	11	Bare Unit / Width (cm)	7
Packed unit / Height (cm)	5	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	16	Bare Unit / Length (cm)	13
Packed unit / Weight	0.2	Bare Unit / Weight (kg)	0.13
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	55
Minimum humidity	5	Maximum humidity	95
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 21:39:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)