

## AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR DISPLAY OF 6 CORAL



Ref. SPG345P



30/08/2024 21:44:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR DISPLAY OF 6 CORAL



Ref. **SPG345P**

**SPENGLER precision in electronic version**

### Description

Display of 6 Autotensio coral wrist.

MEASUREMENT  
ACCURACY• SPENGLER MEASUREMENT•  
WHO risk indicator

EASY &  
CONVENIENT• Fully automatic• 2

buttons• LCD display  
• 120 time-stamped memories•  
Fast and quiet

ADDITIONAL INFORMATION•  
Wrist circumference 13.5-21 cm  
• Comes with 2 LR03 AAA batteries and a storage case

### Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446041601	HS Code	90189010
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.62
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 21:44:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)