

DISPLAY 6 UNITS AUTOTENSIO ARM- ARM BPM ADULT MEDIUM (22X36CM) CORAL



Ref. SPG445P



CE
0459

DISPOSITIF MÉDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2
2 ANS
GUARANTEE

 30/08/2024 21:45:10 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

DISPLAY 6 UNITS AUTOTENSIO ARM- ARM BPM ADULT MEDIUM (22X36CM) CORAL



Ref. **SPG445P**

Ease of use in everyday life for self-measurement family use

Description

Display of 6 Autotensio arms Coral.

HIGH MEASUREMENT
ACCURACY•

SPENGLER MEASUREMENT• Arrhythmia detector•

WHO tricolor risk indicator•

3 average calculations available to monitor the average of morning values EASY TO USE•

Fully automatic• 2
buttons•

New cuff 22-42 cm size M/L

FAMILY

USE• 2 x 60 time-stamped memories

ADDITIONAL INFORMATION

- Comes with 4 LR6 AA batteries and a carrying case• Optional AC adapter
- 2 cuff sizes (S, M/L) available as an option

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446041588	HS Code	90189010
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	3.73
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 21:45:10 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr