

## AUTOTENSIO UPPER ARM (M/L) BLOOD PRESSURE MONITOR SPENGLER CORAL



Ref. SPG445



30/08/2024 21:43:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## AUTOTENSIO UPPER ARM (M/L) BLOOD PRESSURE MONITOR SPENGLER CORAL



Ref. **SPG445**

**Ease of use in everyday life for self-measurement family use**

### Description

HIGH MEASUREMENT  
ACCURACY  
SPENGLER MEASUREMENT • Arrhythmia detector  
WHO tricolor risk indicator  
3 average calculations available to monitor the average of morning values EASY TO USE

Fully automatic • 2 buttons  
New cuff 22-42 cm size M/L

FAMILY USE • 2 x 60 time-stamped memories

ADDITIONAL INFORMATION  
• Comes with 4 LR6 AA batteries and a carrying case  
• Optional AC adapter  
• 2 cuff sizes (S, M/L) available as an option

### Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446041571	HS Code	90189010
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	10
Packed unit / Length (cm)	17	Packed unit / Weight	0.594
Bare Unit / Weight (kg)	0.594		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-25	Max temperature	70
Minimum humidity	5	Maximum humidity	93
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 21:43:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)