

LIAN CUFF/1 TUBE/COTTON/CARTOON DESIGN/INFANT (XS)




Ref. 333039



 30/08/2024 17:47:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LIAN CUFF/1 TUBE/COTTON/CARTOON DESIGN/INFANT (XS)



Ref. 333039

Description

Nylon cuff without coupling for Lian® with 60 cm single tubing bladder for Lian® NM, Lian® Metal, Lian® Nano, Lian® Scenic, Ergo Sphyg®, Lian® Color.
Washable at 30-40° C.

Features

| | | | |
|------------------------------------|---------------|---------------------------|------------|
| CE Marking | NONE | Medical device | Class I |
| Bar code | 3700446000103 | HS Code | 9018901000 |
| Warranty (months) | 6 | Country | Indonesia |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 0.1 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

30/08/2024 17:47:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr