

LIAN CUFF/1 TUBE/COTTON/CARTOON DESIGN/NEWBORN (XXS)




Ref. 333049



 30/08/2024 17:48:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LIAN CUFF/1 TUBE/COTTON/CARTOON DESIGN/NEWBORN (XXS)



Ref. 333049

Description

Nylon cuff without coupling for Lian® with 60 cm single tubing bladder for Lian® NM, Lian® Metal, Lian® Nano, Lian® Scenic, Ergo Sphyg®, Lian® Color.
Washable at 30-40° C.

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446000189	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.055
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 17:48:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr