

SUCCION UNIT ASKIR36, BR LI-ON, 2000ML AUTOCLAVABLE POLYCARBONATE JAR



Ref. BCASPI175



30/08/2024 14:53:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SUCCION UNIT ASKIR36, BR LI-ON, 2000ML AUTOCLAVABLE POLYCARBONATE JAR



Ref. BCASPI175

Designed for professional use and equipped with a powerful suction system that can drain up to 36 l/min max, New Askir 36 LI-ION is lightweight and easy to carry.

Description

Main applications: Ambulance car, Emergency Department, Home care, patient transfer, Room, Tracheostomy, Veterinary Med.

- Proximity function: activation without touch contact

- Three different operating options: The AC/DC adapter, the rechargeable battery and the 12V unit make this surgical Suction Unit a versatile unit guaranteeing CONTINUOUS interventions without the risk of overheating.

- Equipped with a visual and audible alarm system indicating if the battery level is low, the main system is also equipped with a pressure gauge (bar & kPa) and a vacuum regulator.

- Foot with wheels and optional carrying case

- Accessories provided: Makrolon 2858 autoclave polycarbonate fluid collection jar, equipped with an anti-reflux valve with different options / 8 x 14 Silicone hoses (autoclavable) and conical connectors o 10-11-12 mm / Antibacterial & Hydrophobic filter (single patient) / Cannula CH20

- Power cord with Schuko plug, 12V automotive adapter, universal AC/DC

adapter- Engine: Oil-free piston pump requiring no maintenance- Power supply: 14V

with Universal AC / DC adapter / rechargeable internal battery (14.8V 5.2A) Autonomy 70min – Charging time 240min / Automotive adapter (12V 4A)

- Maximum vacuum (adjustable): 0.80 Bar -80 kPa -600 mmHg

- Max free air flow: 36 l/min (without jar)

- Noise level: 60 dBA

- Weight: 3 kg

- Dimensions: 35 x 21 x 18 cm

 30/08/2024 14:53:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SUCCION UNIT ASKIR36, BR LI-ON, 2000ML AUTOCLAVABLE POLYCARBONATE JAR



Ref. [BCASPI175](#)

Features

| | | | |
|------------------------------------|---------------|---------------------------|-----------|
| CE Marking | CE0123 | Medical device | Class IIa |
| Bar code | 8054619101391 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Italy |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 26 | Packed unit / Height (cm) | 30 |
| Packed unit / Length (cm) | 39 | Packed unit / Weight | 5 |
| Bare Unit / Weight (kg) | 3 | | |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | -20 | Max temperature | 45 |
| Minimum humidity | 1 | Maximum humidity | 85 |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

30/08/2024 14:53:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr