

SUCCION UNIT ASKIR36, BR LI-ON, 1000ML FLOVAC DISPOSABLE LINER



Ref. BCASPI176



30/08/2024 14:54:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SUCCION UNIT ASKIR36, BR LI-ON, 1000ML FLOVAC DISPOSABLE LINER



Ref. BCASPI176

Designed for professional use and equipped with a powerful suction system that can drain up to 36 l/min max, New Askir 36 LI-ION is lightweight and easy to carry.

Description

Mobile surgical Suction Unit designed to drain body fluids and perform oral, nasal and tracheal aspiration in adults and children.

Main applications: Ambulance car, Emergency Department, Home care, patient transfer, Room, Tracheostomy, Veterinary Med.

-Proximity function: activation without touch contact- Composed of:

New Askir 36 LI-ON Suction Unit + Flovac 1L kit (1x Flovac 1L reusable jar / 1x Sterile UU pocket 1L Flovac® / 1x Set tubes O8X14mm + connectors / 1x Spread)

- Three different operating options: The AC / DC adapter, the rechargeable battery and the 12V unit make this surgical Suction Unit a versatile unit guaranteeing CONTINUOUS interventions without risk of overheating.

- Equipped with a visual and audible alarm system indicating if the battery level is low, the main system is also equipped with a pressure gauge (bar & kPa) and a vacuum regulator.

- Foot with wheels and optional carrying case

- Accessories provided: Fluid collection jar, equipped with an anti-reflux valve with different options / 8 x 14 Silicone hoses (autoclavable) and conical connectors o 10-11-12 mm / Antibacterial & Hydrophobic filter (single patient) / Cannula CH20

- Power cord with Schuko plug, 12V automotive adapter, universal AC / DC adapter

- Engine: Oil-free piston pump requiring no maintenance- Power supply: 14V

with Universal AC / DC adapter / rechargeable internal battery (14.8V 5.2A) Autonomy 70min – Charging time 240min / Automotive Adapter (12V 4A)

- Maximum vacuum (adjustable): 0.80 Bar -80 kPa -600 mmHg

- Max free air flow: 36 l/min (without jar)

- Noise level: 60 dBA

- Weight: 3 kg

- Dimensions: 35 x 21 x 18 cm

 30/08/2024 14:54:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SUCCION UNIT ASKIR36, BR LI-ON, 1000ML FLOVAC DISPOSABLE LINER



Ref. BCASPI176

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661041008744	HS Code	90189084
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	3
Bare Unit / Weight (kg)	3		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	1	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 14:54:08 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr