

SUCTION UNIT CA-MI NEW HOSPIVAC 350, AUTOCLAVABLE JARS, 2 X 5000 ML- FULL 5 VERSION



Ref. BCASPI506



 30/08/2024 15:13:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SUCTION UNIT CA-MI NEW HOSPIVAC 350, AUTOCLAVABLE JARS, 2 X 5000 ML- FULL 5 VERSION



Ref. BCASPI506

Designed to drain the patient's body fluids, gases, tissues and bones, during or after a procedure. The state-of-the-art oil-free pump system guarantees excellent suction capabilities and maximum vacuum pressure in seconds, without any need for maintenance.

Description

-with foot control and automatic suction inverter

Designed to drain the patient's body fluids, gases, tissues and bones, during or after a procedure. The state-of-the-art oil-free pump system guarantees excellent suction capabilities and maximum vacuum pressure in seconds, without any need for maintenance.

- The precise dashboard and the full range of accessories that come with it make it the ideal device for surgical aspiration and liposuction.
- Accessories included:- Autoclavable

jars with overflow valve

-8x14 Autoclavable Silicone Tubes + Conical connector o 10-11-12mm

-Antibacterial and hydrophobic filter

-Safety container 220 ml

- Technical characteristics

:- Engine: Oil-free piston pump requiring no maintenance

-

Power supply: 230V-50Hz

-

Maximum vacuum (adjustable): -0.90 Bar -90 kPa -675 mmHg

- Max free air flow: 60 l/min

- Noise level: 52 dBA

- Weight: 13 kg

- Dimensions: 46 x 85 x 42 cm

 30/08/2024 15:13:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SUCTION UNIT CA-MI NEW HOSPIVAC 350, AUTOCLAVABLE JARS, 2 X 5000 ML- FULL 5 VERSION



Ref. [BCASPI506](#)

Features

| | | | |
|------------------------------------|---------------|---------------------------|-----------|
| CE Marking | CE0123 | Medical device | Class IIa |
| Bar code | 3661041009932 | HS Code | 90189084 |
| Warranty (months) | 36 | Country | Italy |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 13 |
| Bare Unit / Weight (kg) | 13 | | |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

30/08/2024 15:13:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr