


FLEXIBLE FORCEPS LG : 13CM


Ref. DFPINC01



 30/08/2024 15:29:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

FLEXIBLE FORCEPS LG : 13CM



Ref. DFPINC01

Description

- Crocodile clip
- Flexible
- 13 cm long
- Stainless steel- Ideal for use with a metal speculum (ref. DFSPEC01)
- Intended for veterinary use

Features

CE Marking	NONE	Medical device	No
Bar code	3661040218021	HS Code	90189042
Warranty (months)	12		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	2
Bare Unit / Height (cm)	13	Packed unit / Length (cm)	24
Bare Unit / Length (cm)	13	Packed unit / Weight	0.1
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 15:29:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr