

URINE STRIPS URITOP 11 PARAMETERS REF (100 PCS)

URITOP+®

Ref. DPTDR49

CE



30/08/2024 15:24:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

URINE STRIPS URITOP 11 PARAMETERS REF (100 PCS)

URITOP+®

Ref. DPTDR49

Practical and Fast: 11 Parameters in 60 sec.

Description

- Clinical orientation strips for the detection of 11 parameters in urine
- Multi-parameters for general screening for diabetes, kidney function, or urinary tract infections
- 11 Parameters: GLUCOSE / PROTEINS / pH / BLOOD / KETONE / NITRITES / DENSITY / LEUKOCYTES / BILIRUBIN / UROBILINOGEN / ASCORBIC ACID
- Bottle of 100 strips with scale reading grid
- Results in 60 seconds (except leukocytes: read 90 to 120 seconds)
- Clinical orientation test for professional use or at home
- Optional: Mini URITOP® READER with quantitative and qualitative reading; and storage of results

Features

CE Marking	CE	Medical device	DMDIV (diagnostic In Vitro)
Bar code	3700531701205	HS Code	38220010
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	5	Packed unit / Height (cm)	14
Packed unit / Length (cm)	5	Packed unit / Weight	0.12
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 15:24:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr