

SET OF 20 SUNTECH MONOBLOC ARMBANDS DOUBLE TUBING SINGLE PATIENT SINGLE ADULT LONG (23- 33CM) - SCREW CONNECTOR

SUNTECH

Ref. DSU231086

CE



 30/08/2024 09:37:51 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SET OF 20 SUNTECH MONOBLOC ARMBANDS DOUBLE TUBING SINGLE PATIENT SINGLE ADULT LONG (23- 33CM) - SCREW CONNECTOR



Ref. [DSU231086](#)

Description

Set of 20 Suntech "SOFT" single-patient armbands, royal blue.
Double tubing, adult size 23-33 cm long, screw connection.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446014537	HS Code	9018905040
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.91
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 09:37:51 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr