

MODULE PNI FOR STATION SUNTECH CT40 (ACCESSORY INCLUDED)

SUNTECH

Ref. DSU551400



CE
0413

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
3
3 ANS
GUARANTEE

📅 30/08/2024 08:50:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

MODULE PNI FOR STATION SUNTECH CT40 (ACCESSORY INCLUDED)



Ref. DSU551400

Scalable station ideal for spot measurements of vital parameters

Description

Using the SunTech technology for measuring blood pressure, the Pro VS Check® is the only scalable diagnostic station of this type. Designed to be easily used and adapted in all health care units, the modular design of this device allows you to add blood saturation and thermometry modules tailored to the needs of your patients.

RELIABILITY The manual mode of measuring blood pressure allows reading and checking measurement in real-time, without removing the cuff.

PRECISION The customizable average measurement mode of blood pressure ensures accurate and consistent readings.

CONNECTIVITY Data transfer using IHE protocols with network connection via an integrated Ethernet port or optional WIFI.

FLEXIBILITY The modular design allows you to add optional temperature and SpO2 modules.

SIMPLICITY Intuitive user interface with LCD display.

Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	3700446028800	HS Code	9018905040
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	26	Packed unit / Height (cm)	27
Packed unit / Length (cm)	31	Packed unit / Weight	3.2
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 08:50:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr