

AXILLARY ORAL TEMPERATURE PROBE FOR COVIDIEN FILAC3000 TEMP MODULE

SUNTECH

Ref. DSU551463

CE
0413



 30/08/2024 09:33:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

AXILLARY ORAL TEMPERATURE PROBE FOR COVIDIEN FILAC3000 TEMP MODULE

SUNTECH

Ref. [DSU551463](#)

Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	0120884521008134	HS Code	9018.19.5500
Warranty (months)	3	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	8	Packed unit / Height (cm)	5
Packed unit / Length (cm)	15	Packed unit / Weight	0.086
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 09:33:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr