

**ADAPTOREU 100-240V~/50HZ/28MA AB:  
01/12 OUTPUT: 12V=/300MA/3,6VA**

**KERN®**


Ref. EKB02



 30/08/2024 08:28:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## ADAPTOREU 100-240V~/50HZ/28MA AB: 01/12 OUTPUT: 12V=/300MA/3,6VA



Ref. EKB02

### Description

AC transformer for EKB01

### Features

CE Marking	NONE	Medical device	No
Bar code	3661040201009	HS Code	84231090
Country	Germany		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 08:28:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)