

## PAIR OF ELECTRODS FOR CHILDREN FOR DEFIB. FRED PA-1 SCHILLER



Ref. SSACCE007



30/08/2024 15:44:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## PAIR OF ELECTRODS FOR CHILDREN FOR DEFIB. FRED PA-1 SCHILLER



Ref. **SSACCE007**

### Description

The electrodes of the FRED PA-1 are compatible with all recent models of AED or SCHILLER defibrillator monitors (FRED easyport, FRED® easport plus, DEFIGARD Touch 7, DEFIGARD HD-7), which allows if necessary continuity of care of equipped emergency doctors while keeping the same electrodes.

The positioning of the electrodes on the patient is the same for adults and for children, namely in the anterolateral position.

- Use: children (patients under 8 years of age or weighing less than 25 kg)
- Lifespan: 34 to 36 months
- Cable length: 1.60 m
- Size of the electrodes: 42cm<sup>2</sup> of active surface each

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661041009192	HS Code	90189010
Warranty (months)	36	Country	Austria
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	4
Packed unit / Length (cm)	17	Packed unit / Weight	0.095
Bare Unit / Weight (kg)	0.095		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 15:44:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)