


PAIR OF TRAINING ELECTRODS FOR DEFIB. FRED PA-1 SCHILLER



Ref. SSACCE008



30/08/2024 15:45:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

PAIR OF TRAINING ELECTRODS FOR DEFIB. FRED PA-1 SCHILLER



Ref. **SSACCE008**

Description

- Compatible with FRED PA-1 Trainer and Easy Life training defibrillators used for first aid training
- Reusable and repositionable electrodes
- Designed only for training in the use of a defibrillator in the context of cardiac resuscitation on adults. Do not use on a human.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661041009208	HS Code	90189010
Warranty (months)	36	Country	Austria
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	2
Packed unit / Length (cm)	19	Packed unit / Weight	0.02
Bare Unit / Weight (kg)	0.02		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 15:45:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr