


MANNEQUIN BRAYDEN WITH RED LIGHTS (35X25X60CM)




Ref. SSACCE011



30/08/2024 15:46:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

MANNEQUIN BRAYDEN WITH RED LIGHTS (35X25X60CM)



Ref. **SSACCE011**

Intuitive, easy to clean and long-lasting batteries

Description

This training dummy is equipped with a blood flow display function to help understand and memorize the purpose of cardiopulmonary resuscitation.

Lights indicating chest compression, blood circulation and quality of cardiac massage- Controls the effectiveness of the massage: correct position of the hands, frequency (100 and 120 bpm / min), depth (5 to 6 cm), compression phase and active decompression
 - Comes with: 2 facial skins, 1 artificial lung, 10 insufflation bags, 1 padded carrying bag, 6 LR14 Type C 1.5V batteries.

Features

CE Marking	NONE	Medical device	No
Bar code	3661041009239	HS Code	90230080
Warranty (months)	12	Country	United Kingdom
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	35	Packed unit / Height (cm)	25
Packed unit / Length (cm)	60	Packed unit / Weight	5.42
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 15:46:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
 13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr