


## AC/DC ADAPTOR FOR AIVIA




Ref. SSACCE013



 30/08/2024 15:48:23 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## AC/DC ADAPTOR FOR AIVIA



Ref. **SSACCE013**

### Description

Transformer 24 V DC 75W for cabinet AIVIA

### Features

|                                    |               |                           |           |
|------------------------------------|---------------|---------------------------|-----------|
| CE Marking                         | NONE          | Medical device            | No        |
| Bar code                           | 3661041009253 | HS Code                   | 901889010 |
| Warranty (months)                  | 3             | Country                   | China     |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |           |
| Packed unit / Width (cm)           | 7             | Packed unit / Height (cm) | 6         |
| Packed unit / Length (cm)          | 10            | Packed unit / Weight      | 0.265     |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |           |
| Minimum temperature                | NA            | Max temperature           | NA        |
| Minimum humidity                   | NA            | Maximum humidity          | NA        |
| <b>OTHER FEATURES</b>              |               |                           |           |
| Packed unit / Quantity             | 1             |                           |           |

30/08/2024 15:48:23 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)