

PAIR OF ELECTRODS FOR ADULTS FOR DEFIB. FRED PA-1 SCHILLER



Ref. SSACCE017



30/08/2024 15:52:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

PAIR OF ELECTRODS FOR ADULTS FOR DEFIB. FRED PA-1 SCHILLER



Ref. **SSACCE017**

Description

The electrodes of the FRED PA-1 are compatible with all recent models of AED or SCHILLER defibrillator monitors (FRED easypart, FRED® easport plus, DEFIGARD Touch 7, DEFIGARD HD-7), which allows if necessary continuity of care of equipped emergency doctors while keeping the same electrodes.

The positioning of the electrodes on the patient is the same for adults and for children, namely in the anterolateral position.

- Use: adults (patients aged at least 8 years or weighing more than 25 kg)
- Lifespan: 36 months
- Cable length: 2 m
- Electrode size: 80cm² of active surface each

Features

CE Marking	CE	Medical device	Class I
Bar code	3661041009260	HS Code	90189010
Warranty (months)	36	Country	Austria
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	4	Packed unit / Height (cm)	13
Packed unit / Length (cm)	17	Packed unit / Weight	0.095
Bare Unit / Weight (kg)	0.095		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	1	Max temperature	50
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 15:52:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr