

TRANSPORT BAG FOR FRED PA-1 SCHILLER-STANDARD QUALITY




Ref. SSACCE019



 30/08/2024 15:52:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

TRANSPORT BAG FOR FRED PA-1 SCHILLER-STANDARD QUALITY



Ref. SSACCE019

For safe transport

Description

The standard bag is used to protect and transport the FRED PA-1 or FRED PA-1 Trainer Automated External Defibrillator (AED).

With 2 zippers on each side, it is not necessary to remove the bag to access the defibrillator and the pair of electrodes in case of use.

May contain a defibrillator and a pair of electrodes. The back pocket allows the storage of a pair of additional electrodes and / or consumables.

Semi-rigid cover.

Features

CE Marking	NONE	Medical device	No
Bar code	3661041009284	HS Code	90189010
Warranty (months)	3	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	27	Packed unit / Height (cm)	11
Packed unit / Length (cm)	27	Packed unit / Weight	0.5
Bare Unit / Weight (kg)	0.5		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 15:52:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr