

## LIAN CUFF/1 TUBE/COTTON/CARTOON DESIGN/CHILD (S)




Ref. 333029



 30/08/2024 17:47:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## LIAN CUFF/1 TUBE/COTTON/CARTOON DESIGN/CHILD (S)



Ref. 333029

### Description

Nylon cuff without coupling for Lian® with 60 cm single tubing bladder for Lian® NM, Lian® Metal, Lian® Nano, Lian® Scenic, Ergo Sphyg®, Lian® Color.  
Washable at 30-40° C.

### Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446006877	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.081
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 17:47:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)