

## CLINIC CUFF/1 TUBE/WITH COUPLING # 316 172 AND CLAMP/LARGE ADULT (XL)




Ref. 337870



 30/08/2024 18:05:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## CLINIC CUFF/1 TUBE/WITH COUPLING # 316 172 AND CLAMP/LARGE ADULT (XL)



Ref. 337870

### Description

Single tube Clinic® cuff for Lian® + clamp compatible w/ Welch Allyn, BCI (minition plus) patient monitors and Suntech 247®

### Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028299	HS Code	9018901000
Warranty (months)	6	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.134
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

30/08/2024 18:05:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)