

CLINIC CUFF/1 TUBE/WITH COUPLING # 316 172 AND CLAMP/ADULT (L)




Ref. 337872



 30/08/2024 18:06:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

CLINIC CUFF/1 TUBE/WITH COUPLING # 316 172 AND CLAMP/ADULT (L)



Ref. 337872

Description

Single tube Clinic® cuff for Lian® + clamp compatible w/ Welch Allyn, BCI (minition plus) patient monitors and Suntech 247®

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028282	HS Code	9018901000
Warranty (months)	6	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.095
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:06:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr