

CLINIC CUFF/2 TUBES/WITH COUPLING# 316 151 + CLAMPS/ADULT (M)



Ref. 337951

CE



 30/08/2024 18:07:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

CLINIC CUFF/2 TUBES/WITH COUPLING# 316 151 + CLAMPS/ADULT (M)



Ref. 337951

Description

Dual tube Clinic® cuff + clamp compatible w/ Spengler giant dial Maxi+2® and mercury sphygmomanometers. Marquette Hellige patient monitors

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446000066	HS Code	9018901000
Warranty (months)	6	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.13
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:07:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr