

Ref. 338010



 30/08/2024 18:09:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr



Ref. 338010

Description

Made of yellow Polymex®, with built-in obturating system to allow decontamination by immersion, and Velcro fastening. Legal markings engraved and non-removable.

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446028572	HS Code	9018901000
Warranty (months)	6	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	1
Packed unit / Length (cm)	20	Packed unit / Weight	0.08
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:09:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr