

PAIR OF PINK EARTIPS FOR PULSE II




Ref. 506904

 30/08/2024 18:13:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

PAIR OF PINK EARTIPS FOR PULSE II



Ref. 506904

Description

Pair of ear tips for SPENGLER PULSE II. Standard size, pink color.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030025	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.004
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:13:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr