

PAIR OF BLACK EARTIPS SILICONE (SMALL) FOR CARDIOPRESTIGE II



Ref. 506906



30/08/2024 18:14:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

PAIR OF BLACK EARTIPS SILICONE (SMALL) FOR CARDIOPRESTIGE II



Ref. 506906

Description

Pair of SPENGLER CARDIO PRESTIGE II stethoscope ear tips.
Black color, small size, nano silicone.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446034382	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.004
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:14:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr