

PAIR OF GREEN EARTIPS FOR PULSE II PEDIATRIC



Ref. 506907

 30/08/2024 18:14:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

PAIR OF GREEN EARTIPS FOR PULSE II PEDIATRIC



Ref. 506907

Description

Pair of SPENGLER pediatric PULSE II stethoscope ear tips.
Green color, standard size.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030681	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.004
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:14:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr